

Pipette Tip Packaging Options

Pipette Tip Packaging Options:

Choose the right tip and the right packaging option. Each packaging option is identified in the tip ordering information tables by icons as presented below:

Bulk Tips



Racked Tips



e•dek



TipStation



Platinum Reload System



QuickRack



Gel Loading Tips



Bulk Tips

1,000 loose tips packaged in resealable bags

- Available Tips: Wide selection



Racked Tips

Sorenson tips are available in sturdy racks that withstand several autoclave cycles.

CellRacks hold tips steadily in fully enclosed cells, minimize static and are excellent platforms for single and multichannel pipetting. **Hinged** racks hold the extremely popular 1000XT (50-1250 µl) tips, and traditional 1000 µl standard tips are placed in 100-hole **standard** racks.

- Available Tips: Wide selection



e•dek – Paperboard Rack System

- Easy to recycle – compatible with most paperboard recycling programs
- Used racks require minimal storage space – just flatten and recycle
- Manufactured from paperboard sourced from sustainable forests
- Printed with soy and vegetable based inks
- Less weight equals reduced energy consumption and cost for shipping
- Available Tips: Sterile OneTouch pipette tip range



Pipette Tip Packaging Options



TipStation – Multichannel Pipetting Platform

- Saves space: 192 tips per deck, 5 decks per stack (960 tips)
- Decks designed in 2 x 96 format to match microtiter plate formats
- Stable platform for loading tips
- Can be used with 8 or 12 channel pipettors
- Lid can be positioned to expose only 96 tips
- Available Tips: 10 μ l LongReach and 200 μ l Standardization



Platinum Reload System

- Saves space
- Less waste – environmentally friendly
- Durable, reusable, open-style hinged racks
- No special transfer device needed
- Fit other similar reload system racks
- Clear reload packaging = easy tip identification
- Available Tips: 10 μ l LongReach; 200 μ l Standardization; 300 μ l NX; and 1000 μ l XT (50-1250 μ l)



QuickRack Environmental Packaging

- Quickly loads 96 tips into empty racks
- Compatible with most manufacturers' racks
- Reduces waste
- Available Tips: 10 μ l LongReach; 200 μ l Standardization; 300 μ l NX; and 1000 μ l XT (50-1250 μ l).



Gel Loading Tip Packaging

- Gel loading tip racks (each rack is a single 'pack' of 200 tips) are color coded by outside tip diameter.
 - White racks – Round (0.57 mm O.D.) tips
 - Yellow racks – Flat (0.37 mm O.D.) tips
 - Blue racks – Flat (0.17 mm O.D.) tips
- MultiFlex (1-200 μ l) gel loading tips are available bulk (round) and racked (round and flat)
- MultiFlex E (0.5-10 μ l) gel loading tips are available racked (round and flat)
- MiniFlex (0.1-10 μ l) gel loading tips are available racked (round and flat)